

## KRYOKONTUR KK2+

In the KRYOKONTUR KK2+ version the device is supplied with a mobile swivel arm and shelf



## KRYOKONTUR KK 2

Cryolipolysis with the most innovative technology – MADE IN GERMANY by Ergomed, with more than 15 years' experience in the development and manufacture of cryotherapy equipment.

### TECHNICAL DATA

#### FUNCTIONAL PRINCIPLE

Cooling compressor, CFC and HFC free (especially environmentally friendly)  
Microprocessor control

#### TEMPERATURE RANGE

Optimised temperature profile for effective and safe treatment

#### SERVICING

Annually per item with device check

#### WARRANTY

24 months



#### WEIGHT

approx. 50 kg

#### DIMENSIONS

Depth	710 mm
Height	980 mm
Width	440 mm

#### POWER CONSUMPTION

max. 0,9 kW, 230 Volt, 50 Hz

EC-Certificate acc. 93/42/EEC Annex V • CE 0044 TÜV NORD CERT GmbH • Subject to technical modifications.

The procedure and the device are protected by utility patent. As at 3/2017

[WWW.KRYOKONTUR.DE](http://WWW.KRYOKONTUR.DE)

Manufacturer

INTECCS LTD ERGOMED medical · 44319 Dortmund · Germany

[kontakt@inteccs.de](mailto:kontakt@inteccs.de)

# COOLING TO REMOVE FAT CELLS

## CRYOLIPOLYSIS WITH THE MOST UP-TO-DATE TECHNOLOGY



**KRYO  
KONTUR**

# ABOUT KRYOKONTUR

The KRYOKONTUR system offers completely new options in the treatment of unwanted fat deposits.

It has long been known that under a temperature of 7°C fat cells are damaged.

This process is called cryolipolysis. During this process fat cells are destroyed then broken down and eliminated by a natural metabolic process.

The optimised temperature profile makes treatments with the KRYOKONTUR almost pain-free.

A cooling applicator and a transfer pad in combination with a constant vacuum ensure an even transfer of cold to the surface of the skin. Using this method, depending on the area to be treated, the fat layer is treated for 35-70 minutes. Fat cells are permanently damaged. The metabolisation, decomposition and elimination of the destroyed fat cells, lasts up to 4 months. Apart from a short-term reddening of the skin and slight sensitivity to touch, side effects are so far unknown.

With the KRYOKONTUR technology a device has been developed that meets the high requirements for functionality and simple operation.

KRYOKONTUR is the first device in the world to be operated in an environmentally friendly way without CFCs and HFCs with solid German compressor technology.

Different applicators for the areas of stomach, hips, saddlebags and the insides of thighs ensure optimum treatment options.



DURING TREATMENT



AFTER TREATMENT  
WITH THE 340 CM2 APPLICATOR



TREATMENT OF  
INSIDES OF THIGHS



TREATMENT OF  
SADDLEBAGS AND CELLULITE



## KRYOKONTUR KK2

The KRYOKONTUR KK2 is the first off-the-shelf cryolipolysis device with the quality seal Designed and Made in Germany.

# CRYOLIPOLYSIS WITH KRYOKONTUR MADE IN GERMANY



## PAIN-FREE

The treatment of unwanted fat cells with cold is almost pain-free. As a rule there are no sensations of pain as the cooling has an anaesthetic effect.



## EFFECTIVE

Innovative cooling compressor technology enables the treatment of large areas on all relevant parts of the body. Even temperature distribution and a special temperature profile enable the achievement of optimum results without skin damage.



## INNOVATIVE

Solid technology from Germany, the most up-to-date microprocessor control and a great variety of treatment applicators guarantee a future-oriented equipment system.



## PROFITABLE

Minimal staff requirement as extremely easy to handle - equipment needs no operating during treatment. Only one transfer pad per treatment necessary. Minimal consumption of energy (approx. 20 cents per hour).

## YOUR GAI

- 340 cm2 applicator
- no additional costs apart from transfer pad
- low maintenance
- 24 months warranty
- minimal staff requirement
- low energy costs

