



ERGOMED
medical

KRYODERM V



The successful device for cryocosmetic

For removing age- and pigmentation-marks, warts, keloids, etc.

Local cryotherapy with state-of-the-art technology - energy-saving, environment-friendly, cost-effective.

For dermatology and cosmetics.

KRYODERM V offers an incomparable range of application in the cosmetic field.

Its main operating field is the 100 % removal of age spots, acne, and keloids as well as most kinds of warts and hemangiomas.

KRYODERM V-treatments are painless and don't leave any scars.

Cryo-Cosmetic with KRYODERM V is a biological treatment and easy and safe in success and application.

This device has been developed according to the latest freezing technology. It meets the latest standards of functionality, safety and is environment-friendly.

It is also the first worldwide device which operates FCHC an FHC free with a solid compressor-technology.

The temperature of -35°C is directly obtainable even after numerous treatments.

Technical Data:

Functioning:

compressor, FCHC and FHC-free
(environmentally very friendly)

Temperature range: -35°C bis $+20^{\circ}\text{C}$

Timer:

preset signal clock, 1 to 99 seconds

Size:

depth	400 mm
height	400 mm
width	420 mm

weight: appr. 19,5 kilograms

Power consumption:

appr. 0,4 kW, 230 Volt, 50 Hz

Maintenance:

KRYODERM V is maintenance-free.

Testing certificates:

This machine has been GS-safety approved by TÜV Rheinland –Germany.

EG-Certifikat G3 06 08 60490 001
CE 0123 IEC 601-1/EN 60601-1-2



Treatment



bevor



after

Development and Production

INTECCS LTD
Abtlg. ERGOMED medical
Webershohl 53
D-44319 Dortmund
Germany

Telefon: + 49 (0) 231 610 39 44
Telefax: + 49 (0) 231 610 48 69
E-Mail: kontakt@ergomed-medical.de
Internet: www.ergomed-medical.de



ERGOMED
medical

Products:

KRYODERM -
KRYOLIND -
KRYOKONTUR
KRYODENT -